



LIGHTSHEER

This is how confidence feels

LightSheer® Laser Hair Removal

A collection of Before & After images
from around the world



We don't just build technology, we build confidence.

What we do changes the way people feel about themselves.
We don't just remove unwanted hair, we remove doubts and worries
and the barrier between people and what they can be.
We create, design and manufacture confidence.

This is how **confidence** feels



The LightSheer® Family:

LightSheer® DESIRE™

The LightSheer DESIRE uses an 805nm or 1060nm diode laser wavelength with a variety of handpieces, utilizing two advanced technologies.

The LightSheer DESIRE is simple to operate, treatments are easy to perform and can be scalable to meet your clinic's needs.



LightSheer® DUET™

LightSheer DUET is a complete solution for laser hair removal clinics utilizing two 805nm handpieces, ChillTip and High-Speed.



Laser Hair Removal 805nm

Lower Abdomen



Baseline



Followup after 6 treatments

LightSheer DUET, High Speed handpiece, 805 nm wavelength Skin type III, Fluence: 8.5-10J/cm², Pulse duration: 30 ms

Courtesy of Daniela Cangelosi, MD

Laser Hair Removal 805nm

Back and Shoulders



Baseline



After 2 treatments

LightSheer DUET, High Speed handpiece Fluence: 10J/cm², Pulse duration range: 30-70 ms (Custom mode), 805 nm wavelength

Courtesy of Mitchel Goldman, M.D.

Laser Hair Removal 805nm

Axilla



Baseline



6 months follow up after 3 treatments

LightSheer DUET, High-Speed handpiece Skin type IV, Brown hair , Fluence: 11 J/cm², Pulse duration: 60-66ms, 805 nm wavelength

Courtesy of Suzanne Kilmer, M.D.

Laser Hair Removal 805nm

Lower Abdomen



Baseline



6 months follow up after 5 treatments

LightSheer DUET, High Speed handpiece, 805 nm wavelength

Courtesy of Girish Munavalli, M.D., M.H.S

Laser Hair Removal 805nm

Chest



Baseline



2 months follow-up after 5 treatments

LightSheer DESIRE, High Speed handpiece, 805nm wavelength, Fluence: 7.5-8 J/cm², Pulse Duration: 100 ms

Courtesy of Shweta Virmani, M.D.

Laser Hair Removal 1060nm

Upper Back and Shoulders



Baseline



6 months follow up after 5 treatments

LightSheer INFINITY, High Speed handpiece, 1060 nm wavelength Skin type IV, Fluence: 6-7.5 J/cm², Pulse Duration: 30 ms

Courtesy of Sushil Tahiliani, M.D.

Laser Hair Removal 1060nm

Axilla



Baseline



6 months follow up after 5 treatments

LightSheer INFINITY, High Speed handpiece, 1060 nm wavelength Skin type IV, Fluence: 12 J/cm², Pulse Duration: 50-100 ms

Courtesy of Sushil Tahiliani, M.D.

Laser Hair Removal 1060nm

Axilla



Baseline



6 months follow up after 5 treatments

LightSheer INFINITY, High Speed handpiece, 1060 nm wavelength Skin type IV, Fluence: 12 J/cm², Pulse Duration: 50-100 ms

Courtesy of Sushil Tahiliani, M.D.

Laser Hair Removal 1060nm

Lower Arms



Baseline



6 months follow up after 5 treatments

LightSheer INFINITY, High Speed handpiece, 1060nm wavelength Skin type V, Fluence: 7.5-12.5 J/cm², Pulse Duration: 30-50 ms

Courtesy of Sushil Tahiliani, M.D.

Wrinkle Reducation 1060nm

Under eye



Baseline



3 months follow-up after 3 treatments

LightSheer INFINITY, ET adaptor, 1060 nm wavelength, Fluence 40-50 J/cm²

Courtesy of Lumenis Ltd.

Pigmented Lesions 805nm

Forehead



Baseline



One month follow-up after 3 treatments

LightSheer INFINITY, ET adaptor, 805 nm wavelength, Fluence: 40-45 J/cm², Pulse duration: 20-22.5 ms

Courtesy of Lumenis Ltd.

Pigmented Lesions 805nm

Forehead



Baseline



BL vs 3MFU after 4 treatments

LightSheer, ET adaptor, 805 nm wavelength, Fluence: 50 J/cm², Pulse duration: 25-30 ms, 2 passes

Courtesy of Lumenis Ltd.

Pigmented Lesions 805nm

Forehead



Baseline



BL vs 1MFU after 5 treatments

LightSheer, ET adaptor, 805 nm wavelength, Fluence: 50 J/cm², Pulse duration: 30 ms, 2 passes

Courtesy of Lumenis Ltd.

Your journey with LightSheer has just begun

We want your business to grow with LightSheer. That's why we have put together a support package that begins with the first conversation and will stay with you as your business begins to see the benefits. We provide a support package which ensures that you have all the training you need and receive the most up to date clinical materials. Our direct to consumer marketing activities support the brand, keeping LightSheer the treatment of choice when it comes to laser hair removal.



Lumenis Practice Locator

Thousands of patients search for treatment providers every month. Help them find your clinic by registering your practice on the Lumenis Practice Locator. We run campaigns through influencers, social media and viral videos to promote this directory.



Lumenis online

Consumers have many questions about aesthetic. Our patient-facing websites will give them the re-assurance and confidence they need to move forward with the treatments.



Team Training

We install the system and train your staff, so that you are able to assess skin types and conditions, creating the perfect treatment plan for each individual patient. Our expert clinical training means they will handle consultations and treatments with confidence.



Your Most Valuable Resource

Your dedicated representative is an email, text, or phone call away. We want you to feel part of the Lumenis family, so please feel free to contact us any time. We'll always do our best to get back to you within 24 hours.



Brand Your Clinic

As a Lumenis client you will have access to the Lumenis Physician Portal. Here, you can find materials such as brochures and posters. You will also have access to a wealth of clinical and product information.



Extensive Service Network

We are committed to providing high quality, comprehensive service for medical equipment repair and maintenance. Lumenis service provides you with peace of mind, so you can focus on patient care, rather than on the equipment.



Lumenis - Five Decades of Innovation

**Over 50,000 Systems Installed
in 6 Continents**





AMERICAS

San Jose, CA, USA

T +1 408 764 3000

F +1 408 764 3999

EMEA

Dreieich Dreieichenhain, Germany

T +49 6103 8335 0

F +49 6103 8335 300

Roma, Italy

T +39 06 90 75 230

F +39 06 90 75 269

Hertfordshire, UK

T +44 20 8736 4110

F +44 20 8736 4119

JAPAN

Tokyo, Japan

T +81 3 4431 8300

F +81 3 4431 8301

ASIA PACIFIC

Beijing, China

T +86 10 5737 6677

F +86 10 5737 6767

Kowloon, Hong Kong

T +852 217 428 00

F +852 272 251 51

Gurgaon, India

T +91 124 422 07 95



PB-2002951 Rev C